

# **XPANSION**<sup>®</sup>

MICRO-AUTOGRAFTING KIT



## **Facilitates Durable, Autologous Closure for Complex Wounds:**

- Provides epidermis and dermis to the wound bed
- Offers a solution for high-risk patients

## **Minimizes Size of Donor Site:**

- Can greatly expand the coverage-to-donor-site ratio
- May facilitate faster closure of donor site
- Reduces size of painful donor site

## **Can Be Utilized Outside the OR:**

- Requires only local anesthetic
- Provides an additional option for complex wound management

**Appropriate for a range of complex wounds including:**

**Burns**

**Traumatic Wounds**

**Chronic Wounds**

**Pressure Ulcers**


**Surgical Wounds**



## Why Should You Use Xpansion?

The Xpansion micro-autografting kit allows clinicians to greatly expand the coverage-to-donor-site ratio of autografts, minimizing the size of the patient’s painful donor site while still facilitating durable, autologous closure. The procedure can be performed outside the OR.

## Expand Coverage up to 25 Times the Size of the Skin Graft



The Xpansion device utilizes Meek’s theory of expansion to maximize autologous donor tissue.

Each segment of skin will allow epidermis to migrate from its edges in all directions. By creating micro-autograft segments, migration can occur in multiples allowing greater and greater expansion of a given size of skin graft.

## Ordering Information

Product Name	Catalog Number	Quantity
Xpansion	86EX-PAND-AC	1/box

To order Xpansion or for more information:

✉ [customerservice@appliedtissue.com](mailto:customerservice@appliedtissue.com) ☎ 781-366-3848

